

# SURGITRON<sup>®</sup> DUAL EMC 90

Advanced 4.0 MHz Radiofrequency Technology

**PRECISION • VERSATILITY**



**ellman<sup>®</sup>**  
Experts in Precision Surgery.

# Surgitron® Dual EMC 90 Energy Source

---

With over **50 years of experience**, over 70 patents and more than 200 journal articles, Ellman is your trusted worldwide partner for surgical products and services.

The Surgitron® Dual EMC 90 unit represents advanced radiofrequency technology that provides unparalleled surgical control, **precision, versatility** and **safety**. The high frequency of 4.0 MHz minimizes heat dissipation and thus cellular alteration.


## Five Distinct Waveforms for Optimum Results

---

### 1. Fully Filtered (Cut)



- Micro-smooth cutting
- Negligible lateral heat
- Minimal cellular destruction
- Best cosmetic results. Fastest healing
- Ideal for skin incision and biopsy
- 4.0 MHz



### 4. Fulguration




- Maximum hemostasis
- Ideal for intentional tissue destruction
- 4.0 MHz



### 2. Fully Rectified (Cut/Coag)




- Cutting with hemostasis
- Ideal for sub-cutaneous tissue dissection and planing. Especially useful in vascular areas where producing minimal amounts of lateral heat and tissue damage
- 4.0 MHz




### 5. Bipolar


- Pinpoint, micro-coagulation
- Minimal charring or tissue necrosis
- Ideal for coagulation in and around critical anatomy
- 1.7 MHz



### 3. Partially Rectified (Hemo)



- Coagulation / Shrinkage
- Ideal hemostasis with controlled penetration
- For cutting with maximum hemostatic control
- 4.0 MHz



### Clinical Citations

1. Olivar, A.C., et al, Ann Clin Lab Sci. (1999); 29(4): p281-285.
2. Data on file.
3. Niamtu, J., Chapter 4B, "Radiowave Surgery in Oral and Maxillofacial Surgery", in Bell, W., *Distraction Osteogenesis of the Facial Skeleton*, 2007, p30-37.
4. Bridenstine, J.B., Derm Surgery (1998); vol 24, p397-400.
5. Ericsson, E., et al, The Laryngoscope (2007); vol 117, p654
6. Merz, Derm Surgery (2002); vol 28, p735-738
7. Eremia, S., et al, Dermatol Surg (2001); 27: p1052-1054.

## Distinct Benefits for Your Practice and Your Patients



- **Precision** – create precise incisions in a variety of tissue structures
- **Versatility** – no other energy-based technology has the surgical versatility of Ellman
- **Quick Recovery** – with less tissue destruction, healing is hastened and your patients can recover quickly
- **Decreased Post-Operative Pain** 4.0 MHz causes less trauma
- **Decreased Post-Surgical Edema** low temperature equals less tissue destruction
- **Less Burning or Charring of Tissue** 4.0 MHz minimizes burning of tissue, unlike laser or conventional lower frequency electrosurgery

## Features

- **Dual Frequency** combining two optimized frequencies Monopolar (4.0 MHz) and Bipolar (1.7 MHz) for maximum surgical benefits – for outstanding precision and control
- **Digital Control Panel** facilitates easy operation and a clear view of settings
- **Solid State Circuitry** for dependable and consistent energy emission
- **Safety Indicators** provide visual and auditory alerts
- **Parameter Recall** allows rapid set-up for subsequent procedures



## 4.0 MHz Minimizes Lateral Thermal Damage and Maximizes Precision



Reference - Olivar, AC, et al, Ann Clin Lab Sci. 1999 Oct-De;29(4):

Source: Golio, JM, et al, "RF and Microwave Applied Systems", The RF and Microwave Handbook, 21-2.

- Ellman radiofrequency technology produces one-third the lateral thermal damage as compared to conventional electrosurgery
- Ellman radiofrequency technology produces one-half to one-third of the lateral thermal damage versus most lasers

## Shown with Surgitron® Dual EMC 90, Surg-e-Vac™ and Cart

### Surgitron® Dual EMC 90 Specifications

<b>Product Code</b> EC3A-S30-90	
<b>Description:</b> Surgitron® 4.0 Dual R790 IEC	
<b>Dimensions</b>	<b>Output frequency</b>
Height: 5 inches	4.0 MHz Monopolar
Width: 9 inches	1.7 MHz Bipolar
Depth: 13 inches	<b>Line Frequency</b>
Weight: 18 lbs	50 - 60 Hz
<b>Line Voltage</b>	
110/120/220/240 volts	
<b>Output Power</b>	
Monopolar Cut: 90 Watts	
Monopolar Cut/Coag: 65 Watts	
Monopolar Hemo: 45 Watts	
Monopolar Fulgurate: 35 Watts	
Bipolar: 90 Watts	



### How Our Patented Radiofrequency Technology Works

#### Cellular Radiofrequency Absorption

